

Probing Parasitic Gas Evolution in Sodium Ion Batteries

Sarat Alabidun¹, Bethan J, V Davies², Maria Crespo-Ribadeynera¹, Ifan E. L Stephens^{2*}, Magda Titirici^{1*}

Imperial College
London

¹Imperial College London Department of Chemical Engineering, Exhibition rd South Kensington. SW7 2AZ


²Imperial College London Department of Chemistry, White City Campus, 10 wood ln White City W12 7TA



Motivation


- Sodium ion batteries are considered a viable alternative to lithium ion batteries
- These batteries degrade by decomposition of the electrolyte and gas evolution hence, studying the gas evolution mechanism is important
- In-situ and In-operando mass spectroscopy techniques are a great way to study this degradation
- However, current techniques cannot detect very small quantities of volatile species and are therefore not very representative of battery degradation mechanism

Solution- Electrochemistry Mass Spectroscopy (EC-MS) on a Chip!



- Developed and used for aqueous electrochemistry, recently adapted for non-aqueous electrochemistry and use in batteries.(1)(DB Thornton et al)
- Picomole sensitive technique, allowing for detection of very small amounts of volatile species previously undetected in batteries and above current state of the art.


Benchmarking EC-MS Cell compared to Coin Cell



Electrolyte- NaPF₆ EC:DEC
Working Electrode- Hard Carbon
Counter Electrode- Na Metal

- Optimised the designed cell for EC-MS to ensure reproducibility compared to the coin cell benchmark for G1500 hard carbon against sodium metal.(2)
- Oxidation and reduction peaks (also considered sodium insertion and extraction peaks) observed for both cells at around 0V.(2)
- Onset of electrolyte reduction peak observed for both cells close to 1V

Gas Evolution In Sodium Ion Battery



- Maximum 10μmol/min C₂H₄ (This work) detected vs 1μmol/min max.(3)
- Maximum 1μmol/min CO detected(This work) vs 0.2μmol/min max.(4)

CO Corresponding to decomposition of ethylene carbonate to form the SEI.




CO₂ formed at open circuit potential while the cell is at rest. No CO₂ observed when cell begins to cycle possibly due to decomposition of CO?



H₂ formation Due to decomposition of trace H₂O in contact with Na metal,
 $\text{H}_2\text{O} + \text{e}^- \longrightarrow \text{OH}^- + \frac{1}{2} \text{H}_2$

C₂H₄ formed at full discharge (0V) attributed to decomposition of ethylene carbonate



Left- Potential, Current and Gas Flux against Time
Right- Current against Potential

Electrolyte- NaPF₆ EC:DEC
Working Electrode- Hard Carbon
Counter Electrode- Na Metal

Conclusion

The EC-MS is an excellent technique to study gas evolution in battery systems, and has proven to be a few orders of magnitude more efficient in detecting volatile species previously undetected in other mass spectrometry systems.

Acknowledgement

We would like to acknowledge the Damian Cummins scholarship, for providing PhD funding, and Spectro-Inlets for their support with the EC-MS

SpectroInlets
Enabling real-time analysis

REFERENCES

- Trimarco DB, Scott SB, Thilsted AH, Pan JY, Pedersen T, Hansen O, et al. Enabling real-time detection of electrochemical desorption phenomena with sub-monolayer sensitivity. *Electrochimica Acta*. 2018;268:520-30.
- Au H, Alptekin H, Jensen ACS, Olsson E, O'Keefe CA, Smith T, et al. A revised mechanistic model for sodium insertion in hard carbons. *Energy & Environmental Science*. 2020;13(10):3469-79.
- Zhang L, Tsolakidou C, Mariyappan S, Tarascon J-M, Trabesinger S. Unraveling gas evolution in sodium batteries by online electrochemical mass spectrometry. *Energy Storage Materials*. 2021;42:12-21.
- Goktas M, Bolli C, Berg EJ, Novak P, Pollok K, Langenhorst F, et al. Graphite as intercalation electrode for sodium-ion batteries: electrode dynamics and the missing solid electrolyte interphase (SEI). *Advanced Energy Materials*. 2018;8(16):1702724